

EXHIBIT A

JUL 22 2004 9:16AM

NO. 346 P. 2/3

Biomechanical Testing of Tibialis Anterior Graft Tibial Tunnel Fixation with Bioabsorbable RetroScrews and BioScrew XtraLok in Porcine Bones

H C Chang, J Nyland, A Nawab, R Borden, D N M Caborn

Division of Sports Medicine,
Department of Orthopaedic Surgery,
University of Louisville,
Louisville, Kentucky 40202, USA.

Haw Chong Chang MBBS, FRCSEd(Orth)
John Nyland, Ed.D., P.T., S.C.S., A.T.C
Akbar Nawab, MD
Robert Borden, MEng, EIT
David N.M. Caborn, MD

Purpose: This study evaluated the failure mode, maximum load at failure, displacement at failure, and stiffness differences of doubled tibialis anterior graft-tibial tunnel fixation using retrograde bioabsorbable interference screws (Arthrex, Naples, FL) and 35-mm BioScrew XtraLok (Linvatec, Largo, FL) after cyclical loading.

Type of Study: Experimental laboratory biomechanical study.

Hypothesis: There is no difference in maximum load, displacement and stiffness at failure of doubled tibialis anterior graft-tibial tunnel fixation using Retrograde (Arthrex, Naples, FL) and 35-mm BioScrew XtraLok (Linvatec, Largo, FL) after cyclical loading.

Methods: Twelve specimens of porcine tibias were divided into 6 matched pairs based on bone mineral densitometry. Wilcoxon tests comparisons were used to assess group differences ($P < .05$).

Results: Maximum load at failure after cyclic loading for the RetroScrew was 778.7 ± 177.5 N, with a displacement of 5.3 ± 2 mm and a stiffness modulus of 204.3 ± 52.9 N/mm. Maximum load at failure after cyclic loading for the BioScrew-XtraLok screw

- 1 -

JUL 22 2004 9:16AM

NO. 346 P. 3/3

was 1436.3 ± 331.3 N, with a displacement of 5.9 ± 2.6 mm and a stiffness modulus of 323.6 ± 56.8 N/mm. Fixation using XtraLok screws displayed greater maximum load at failure than RetroScrew fixation ($P = .028$) as well as greater stiffness ($P=0.046$). Significant differences were not evident for displacement at final pullout. All constructs failed by graft pullout.

Conclusions: Fixation using a single 35-mm BioScrew XtraLok screw displayed increased maximum load at failure and stiffness compared with the 20-mm RetroScrew with 17-mm cortical backup fixation.